



















## Your Complete **Insulin Pen** Reference Guide From *Diabetes Health*

Manufacturer	Products	Comments	For More Information	
<b>Prefilled Insulin Pens (Disposable)</b>				
  <b>Eli Lilly and Company</b>	<b>Humalog Mix 75/25 Pen, 3 ml, 300 units</b>  <b>Humalog Pen, 3 ml, 300 units</b> 	Dial your dose in 1-unit increments, up to 60 units, by turning the easy-to-read adjustable dose knob. These pens are available by prescription only.	<a href="http://www.lillydiabetes.com">www.lillydiabetes.com</a> (800) 545-5979	
	<b>Humulin 70/30 Pen, 3 ml, 300 units</b>  <b>Humulin N Pen, 3 ml, 300 units</b> 	Dial your dose in 1-unit increments, up to 60 units, by turning the easy-to-read adjustable dose knob.		
	  <b>Novo Nordisk</b>	<b>Novolin R InnoLet, 300 units</b> <b>Novolin N InnoLet, 300 units</b> <b>Novolin 70/30 InnoLet, 300 units</b> 		InnoLet is a disposable doser whose dosing dial has the largest, most legible numbers of any insulin delivery system. The large grip and ergonomic design of InnoLet make it easy and comfortable to hold. InnoLet rests securely against the skin, to help minimize tremor or movement at the needle tip and to prevent you from injecting too deeply. InnoLet is designed for use with NovoFine disposable needles.
		<b>NovoLog Mix 70/30 FlexPen, 3 ml, 300 units</b> <b>NovoLog FlexPen, 3 ml, 300 units</b>  <b>NovoLog, 300 units</b>  <b>Levemir, 300 units</b> 		Both of these disposable pens have large-dose dialing windows and allow dosing in 1-unit increments, up to 60 units. They also allow easy dose correction by dialing backward or forward. The dose scales reset to 0 after injection, ready for the next dose. The FlexPens are designed for use with NovoFine disposable needles. These pens are available by prescription only.
<b>NovoPen Junior</b> 		NovoPen Junior features a dosing window and a dial-dosing mechanism. The user can administer 1 to 35 units of insulin in 1/2-unit increments. Available in two bright colors, NovoPen-Junior is designed for use with Novo Nordisk insulins in PenFill 3-ml cartridges and with NovoFine disposable needles.		




## Your Complete **Insulin Pen** Reference Guide From *Diabetes Health*

Manufacturer	Products	Comments	For More Information
 <b>Novo Nordisk</b>	<b>NovoPen3</b> 	NovoPen3 features a dosing window and a dial dosing mechanism. The user can administer 2 to 70 units of insulin in 1-unit increments. NovoPen3 is designed for use with Novo Nordisk insulins in PenFill 3-ml cartridges and with NovoFine disposable-needles.	<a href="http://www.novonordisk-us.com">www.novonordisk-us.com</a> (800) 727-6500
	<b>NovoPen3 PenMate</b> 	NovoPen3 PenMate is an automatic needle insertion device designed for quick, easy injections. It is available for use with NovoPen3 and as part of the NovoPen Junior Diabetes Management System.	

### Reusable Insulin Pens That Use Cartridges

 <b>Owen Mumford, Inc.</b>	<b>Autopen</b> 	The Autopen is a multidose injection device available in models that use either 1.5-ml or 3-ml insulin cartridges. The Autopen features a dose selector with large, easy-to-read unit numbers and audible clicks; a side-mounted activation button for easier and gentler automatic delivery; and a clear cartridge holder and dose completion indicator. The two Autopen models for 1.5-ml cartridges offer doses either in 1-unit increments, up to 16 units, or in 2-unit increments, up to 32 units. The two Autopen models for 3-ml cartridges offer doses either in 1-unit increments, up to 21 units, or in 2-unit increments, up to 42 units.	<a href="http://www.owenmumford.com">www.owenmumford.com</a> (800) 421-6936
 <b>sanofi aventis</b>	<b>OptiClik® Pen</b> 	OptiClik is a reusable insulin delivery device (insulin pen) for exclusive use with 3-ml Lantus® (insulin glargine [rDNA origin] injection) or Apidra® (insulin glulisine [rDNA origib] injection) 3-mL Cartridge System. OptiClik® is available in two different colors. For patients that use two different types of insulin, Lantus® or Apidra®, it is recommended to use a different color pen for each insulin. Check the label on your Cartridge System to make sure you have the correct insulin before injecting.	<a href="http://www.opticlik.com">www.opticlik.com</a> (800) 633-1610

### Insulin Pen Needles

 <b>BD &amp; Company</b>	<b>BD Ultra-Fine III Mini Pen Needle</b> <b>BD Ultra-Fine III Short Pen Needle</b> <b>BD Ultra-Fine III Original Pen Needle</b>	BD Ultra-Fine Pen Needles fit all diabetes pens and dosers sold in the U.S. BD offers three different needle sizes to meet all patient needs. Some of the shortest and thinnest pen needles in the world. 	<a href="http://www.BDdiabetes.com">www.BDdiabetes.com</a> (888) 232-2737
 <b>Owen Mumford, Inc.</b>	<b>31G x 6 mm ultra short</b> <b>31G x 8 mm short</b> <b>29G x 12 mm original</b>	<ul style="list-style-type: none"> <li>• Three sizes to suit the needs of all patients</li> <li>• Fit quickly and easily to all major insulin delivery devices</li> <li>• FREE pen needle remover in every box to make removing pen needles safe and convenient.</li> </ul>	<a href="http://www.owenmumford.com">www.owenmumford.com</a> (800) 421-6936