

Your Complete **Infusion Sets** Reference Guide From *Diabetes Health*

Manufacturer	Brand Name	Connector Type	Cannula Type	Available Cannula Lengths	Available Tubing Length (in inches)	Available Automatic Inserter?	Disconnect at Site?	For More Information
Fresenius Kabi	ACCU-CHEK Rapid-D	Luer-lock	28-gauge surgical stainless steel needle 90° insertion angle	6 mm 8 mm 10 mm	24, 31, 43	No	QR*	(800) 280-7801 www.disetronic-usa.com www.chooseyourset.com
Dana Diabecare	Dana	Dana proprietary	Straight needle	7 mm 9 mm	27, 43	No	QR*	(866) 342-2322 www.danapumps.com
	Dana	Dana proprietary	Plastic cannula	7 mm 9 mm	27, 43	No	QR*	
Disetronic Sterile Products	ACCU-CHEK Ultraflex	Luer-lock	Soft Teflon cannula 90° insertion angle	8 mm 10 mm	24, 31, 43	Yes	QR*	(800) 280-7801 www.disetronic-usa.com www.chooseyourset.com
FIFTY50	FIFTY 50 THINSet	Luer-lock, Works with all pumps* including Medtronic MiniMed Paradigm	90o Soft Teflon cannula	6mm 9mm	23, 43	Yes	Yes	(800)746-7505 www.FIFTY50.com
Medtronic/ MiniMed	Quick-set	Luer-lock, Paradigm	Soft Teflon	6mm 9mm	23,43	Quick-serter	Yes	(800) 646-4633 www.medtronic.com
	Silhouette	Luer-lock, Paradigm	Soft Teflon	13mm 17mm	23,43	Sil-serter	Yes	
	Sof-Set	Luer-lock, Paradigm	Soft Teflon	6 mm 9 mm	24, 42	Sof-serter	No	
	Sure-T	Luer-lock	Steel Needle	6 mm	23	No	No	
	Polyfin	Luer-lock, Paradigm	Steel Needle	14 mm	24, 42	No	No	
Smiths Medical MD	Cleo 90	Luer-lock	Teflon cannula 90° insertion angle	6 mm 9 mm	24, 31, 42	Integrated insertion/retraction device	Yes	(800) 826-9703 www.cleoinfusionsets.com
SimpleChoice	Twist	Luer-lock	Teflon soft catheter 30° insertion angle	12 mm 17 mm	23, 43	Yes, with adapter	Yes	(877) 773-2879 www.mysimplechoice.com
	Twist PRO	Paradigm	Teflon soft catheter 90° insertion angle	9 mm	23, 43	Yes, with adapter	Yes	
	Easy	Luer-lock	Teflon soft catheter 30° insertion angle	12 mm 17 mm	23, 43	Yes, with adapter	Yes	
	Easy PRO	Paradigm	Teflon soft catheter 90° insertion angle	9 mm	23, 43	Yes, with adapter	Yes	
Unomedical	Comfort Comfort Short	Luer-lock	Soft cannula	17 mm 13 mm	23, 31, 43	No	Yes	www.infusion-set.com www.unomedical.com
	Inset	Luer-lock	Soft cannula	6 mm 9 mm	23, 43	Integrated inserter	Yes	
	Contact Detach	Luer-lock	Stainless steel	6 mm 8 mm 10 mm	23, 31, 43	No	Yes	
	insuflon	Luer-lock	Soft cannula	17 mm	23, 31, 43	N/A	N/A	

*"QR" means "quick-release." The set is disconnected from the tubing a couple of inches from the infusion site, leaving a short "tail" of tubing while disconnected.